

3.2.4.1 Salmonid Ceratomyxosis Methods Appendix 1

Jerri L. Bartholomew

Department of Microbiology
Nash Hall 220, Oregon State University
Corvallis, OR 97331-3804
541-737-1856
Jerri.Bartholomew@oregonstate.edu

Molecular Detection of *Ceratomyxa shasta* Using the Polymerase Chain Reaction Assay (Palenzuela et al. 1999)

A. Sample Collection

Because of the specificity of the polymerase chain reaction (PCR) assay, contamination of the samples by bacteria or other contaminants is not a problem and sterility or aseptic techniques are not strictly necessary. However, the risk of cross contamination and carry-over of DNA from other samples is high and sampling must be done with this in mind.

Samples should be collected as follows:

1. Use disposable material for each fish or alternatively (see notes below) destroy the DNA from the dissecting tools between fish. Process one fish at a time.
2. Remove the intestine (or a piece) from the fish and place it on a clean, disposable surface (a piece of aluminum foil works just fine).
3. Cut a small piece of the gut (25 to 100 mg, about 2 to 5 mm), preferably from the lower tract. We use disposable, nonsterile razor blades to cut this small piece, and a simple toothpick to transfer it to the screw-cap vial with 500 μ L DNA extraction buffer (see next section for recipe). Conventional dissecting tools can be used but only if they are completely decontaminated between samples.
4. Dispose of the working surface, change gloves if they came in contact with the tissues, and proceed to the next fish.

Notes on decontamination of material: Conventional methods for disinfecting dissecting tools, such as spraying with ethanol and flaming, not only **DO NOT** completely destroy the DNA, but

3.2.4.1 Salmonid Ceratomyxosis Methods Appendix 1 - 2

they can actually fix it in the material. Procedures which are used to clean the material between samples include:

5. Use of the commercial products DNA AWAY or RNase AWAY. Both destroy DNA quickly and safely. It is sufficient to wipe the material with a paper towel wetted with the product and then rinse it with fresh, distilled water from a lab bottle. These products are not recommended for use with metallic material, but we have used them liberally on forceps and scissors and have not encountered problems.
6. First wipe off any pieces of tissue, blood, or mucus from the tools, then immerse them in 10% bleach (Clorox) solution for *at least* one minute. Rinse before use with distilled water from a wash bottle. This concentration of bleach is corrosive for metals after long exposures, but this is not a problem if the exposure is limited to the sampling and tools are rinsed well in distilled water before storage. Use two containers (per person) with bleach during sampling, each containing a set of dissecting tools so that tools can soak between uses. After use, the tool is wiped clean with a piece of paper and put back in the bleach solution.

B. Protocol for the Diagnosis of *Ceratomyxa shasta* in Fish Intestines by PCR

1. DNA Extraction

a. Reagents:

i. DNA Extraction Buffer

The buffer is NaCl 100 mM, Tris-HCl 10 mM, EDTA 25 mM, SDS 1%.

To prepare 1 Liter of Extraction Buffer:

into 500 mL HPLC H₂O

20 mL 5M NaCl

5 mL 1M Tris-HCl, pH 7.8

50 mL 0.5M EDTA, pH 8.0

50 mL 20% SDS

pH solution to 8.2 then add HPLC H₂O to a final volume of 1 Liter

Aliquot extraction buffer into 50 mL tubes and seal with parafilm. The aliquots are stable at room temperature for 2+ years.

ii. Proteinase K

Obtain commercially as a stable liquid solution (store at 4°C) or make at a concentration of 20 mg/mL and freeze at -20°C in 1 mL aliquots

iii. RNase A

100mg/mL stock commercially available. Dilute for use as needed.

b. Methods:

When processing tissue, special care must be taken as carry-over of DNA between samples occurs easily. Use new blades and/or tools when handling the tissues. Samples can be fresh or frozen. There is no need to clean the tissues of fat or adherent tissue.

3.2.4.1 Salmonid Ceratomyxosis Methods Appendix 1 - 4

For **10 mL** of master mix:

Sterile molecular grade water	5.9 mL
5x colorless Taq PCR Buffer (Promega)	2 mL
MgCl ₂ 25 mM (Promega)	800 µL
dNTPs (at 100mM; Promega)	80 µL (20 µL each A,T,G,C)
Primer F (CS1 at 100 µM)	50 µL
Primer R (CS3 at 100 µM)	50 µL
Rediload (Invitrogen)	500 µL

Make the mix on ice and under nuclease-free conditions. Aliquot to avoid freeze-thaw cycles and store at -20°C.

3. PCR Protocol

Thaw an aliquot of master mix (18.8 µL per sample) and add 1 U Taq polymerase per sample (usually 0.2 µL, depending on supplier). Assemble the reaction mixture on ice in 20 µL reactions as follows:

- 19 µL master mix
- 1 µL diluted template (sample DNA)

Briefly mix and zip-spin reaction mixture.

PCR program for amplification of *C. shasta*:

95°C/3 min + 34 cycles x {94°C/1 min + 56°C/30 sec + 72°C/1 min} +72°C/10 min.

Analyze 10 µL of each reaction in a 1% agarose gel stained with SYBR Safe DNA stain (Invitrogen) or ethidium bromide. Positive samples will have an amplicon of 640 bp.